

## **Respiratory Viral Panel**

For in vitro diagnostic use

### **CERNER ORDERABLE**

RESPIRATORY VIRAL PCR; RVP PCR 3380

### **CPT CODE**

87633

### **CLINICAL UTILITY**

Respiratory viruses cause acute local and systemic illnesses that range in severity, with the potential to cause severe disease especially in the young and elderly. The NxTAG™ Respiratory Pathogen Panel (RPP) is a qualitative nucleic acid multiplex test intended for the simultaneous detection and identification of 20 viral and bacterial targets in individuals suspected of respiratory tract infections<sup>1</sup>.

### **METHODOLOGY**

Multiplex RT-PCR/Array bead hybridization

### **SPECIMENS**

Nasopharyngeal swabs (NPS) in Universal Transport Medium (UTM)\*

Bronchoalveolar lavage (BAL) collected in a sterile container\*\*

Sputum collected in a sterile container\*\*

### **SPECIMEN STABILITY**

Refrigerated: up to 72 hours

### **SHIPPING**

Ship specimen on ice packs or if greater than 72 hours ship frozen on dry ice

### **CAUSES FOR REJECTION**

Specimens other than those listed above.

### **SPECIFICITY**

Detects 20 Respiratory targets: Influenza A (all subtypes, Influenza A subtype H1, Influenza A subtype H1N1, Influenza A subtype H3), Influenza B, Respiratory Syncytial Virus (RSV) subtype A, RSV subtype B, Coronavirus (229E, OC43, NL63, HKU1), Parainfluenza 1, Parainfluenza 2, Parainfluenza 3, Parainfluenza 4, Human Bocavirus, Human Metapneumovirus, Rhinovirus/Enterovirus, Adenovirus, *Chlamydomphila pneumoniae*, and *Mycoplasma pneumoniae*.

### **ASSAY RANGE**

Qualitative results (Detected/Not Detected) for: Influenza A, Influenza A subtype H1, Influenza A subtype H1N1, Influenza A subtype H3, Influenza B, Respiratory Syncytial Virus (RSV) subtype A, RSV subtype B, Coronavirus 229E, Coronavirus OC43, Coronavirus NL63, Coronavirus HKU1, Parainfluenza 1, Parainfluenza 2, Parainfluenza 3, Parainfluenza 4, Human Bocavirus, Human Metapneumovirus,

1. Reference information can be found in the Indiana University Health Molecular Assay Procedures.

Rhinovirus/Enterovirus, Adenovirus, *Chlamydomonas pneumoniae*, and *Mycoplasma pneumoniae*.

**TURNAROUND TIME**

Daily, 24-72 hours

\*NPS has been cleared by the FDA for use in the RVP assay.

\*\*In-house validation was performed to establish as suitable specimen type

1. Reference information can be found in the Indiana University Health Molecular Assay Procedures.