



Order of Draw

Blood Cultures

1. Aerobic first
2. Anaerobic second



Clear 3ml discard tube

**Use prior to collecting coagulation testing (Blue tube below) & or using butterfly needle.*



Blue (Sodium Citrate)



Gold (SST)

Red

Royal-Blue (Metals)



Orange (RST)



Green (Sodium Heparin)

Light Green (Lithium Heparin PST)



Lavender (EDTA)

Pearl (PST)

Tan (Na-Heparin)

Royal-Blue (Di-Na EDTA)



Grey (Na-Fluoride K+ Oxalate)



Yellow (ACD)