

Order of draw

Blood cultures: 1. Aerobic; 2. Anaerobic



Clear 3ml discard tube

*Use prior to collecting coagulation testing
(Blue tube below) & or using butterfly needle.



Blue (Sodium Citrate)



Gold (SST), Red, Royal-Blue (Metals)



Orange (RST)



Green (Sodium Heparin)



Light Green (Lithium Heparin PST)



Lavendar (EDTA)



Pearl (PST)



Tan (Na-Heparin)



Royal-Blue (Di-Na EDTA)



Grey (Na-Fluoride K+ Oxalate)



Yellow (ACD)

