



**Specimen Transport Temperature Change for Coccidioides
Antigen Detection by EIA, Histoplasma Antigen QN, and
Blastomyces dermatitidis Antigen EIA**

Effective Date:

February 25, 2021

Contact: Sendout Department

Update Description:

As of February 25th, tests for **Coccidioides Antigen Detection by EIA, Histoplasma Antigen QN, and Blastomyces dermatitidis Antigen EIA** to MiraVista must be shipped **FROZEN**. This includes all sample types for each test.

The transport temperature on the label will be updated to “SO FR MISC” and will be updated in the Test Directory.

Please disseminate this information to those who need to know.

Questions:

If there are any questions regarding this test, please contact the IU Health Sendout Department at 317.491.6805. For additional questions or concerns, please contact Kim Kinney, Sendout Department Manager, at 317.491.6850, or email at kjkinney@iuhealth.org