



Laboratory Services

STAT Courier and Time Sensitive Test Menu for Outreach Clients

Below is the list of tests that can be ordered STAT Courier by an IU Health Laboratory Services outreach client. To order a STAT Courier test, make sure to follow the correct STAT test procedure, which includes placing the specimen in a STAT bag.

If there are any questions, please contact Client Services at 317.491.6000 or your account executive.

*For tests marked as “Time Sensitive,” please call Client Services (317.491.6000) for pick-up details. Depending on the time the order is placed (in relation to your scheduled courier pickup time), a STAT Courier may or may not be dispatched.

Test Code	STAT TEST	Testing Department	STAT Courier or Time Sensitive	Time Limit
6003	Ammonia	Chemistry	STAT	N/A
All	Bilirubin	Chemistry	STAT	N/A
All	Blood Bank Testing	Blood Bank	STAT	N/A
All	Blood Gas	Chemistry	STAT	N/A
3132	Blood Parasite Exam	Hematology and Micro	STAT	N/A
5027	Carboxyhemoglobin	Chemistry	STAT	N/A
150	CSF Count	Hematology	STAT	N/A
10029	CSF Cryptococcal Ag	Immunology	STAT	N/A
3024	CSF Gram Stain	Microbiology	STAT	N/A
545	D Dimer	Coag.	STAT	N/A
6015	Electrolytes	Chemistry	STAT	N/A
6716	Ethanol	Chemistry	STAT	N/A
6853	Ethylene Glycol	Chemistry	STAT	N/A
7016	Fetal Fibronectin	Chemistry	STAT	N/A
6025 / 6082	Glucose	Chemistry	STAT	N/A
5	Hematocrit	Hematology	STAT	N/A
4	Hemoglobin	Hematology	STAT	N/A
6155	Ionized Calcium	Chemistry	STAT	N/A
6017	Potassium SerPI QN	Chemistry	STAT	N/A
6846	Volatiles	Chemistry	STAT	N/A
8902	Muscle Biopsy	AP	Time Sensitive	3 Hours
450	Platelet Aggregation <i>(must be scheduled prior)</i>	Special Coag	Time Sensitive	3 Hours
586	Platelet Function Analyzer (PFA-100)	Special Coag	Time Sensitive	4 Hours
587	Platelet Inhibition (P2Y12)	Special Coag	Time Sensitive	3 Hours
6062	Synovial Fluid Cell Count	Hematology	Time Sensitive	4 hours